

Frequently Asked Questions

- **Will ST80i interface with other exercise devices than those listed?**
 - At this time, the ST80i is only compatible with the exercise devices listed. We continue to work on validating additional devices..
- **Will the ST80i interface with all Echo systems for Stress Echo?**
 - Any Echo system that accepts either Analog ECG input or ECG Sync Pulse input will work with ST80i for Stress Echo. The interface cable must be provided by the Echo system manufacturer.
- **Can the ST80i work with rechargeable batteries?**
 - We do not recommend using rechargeable batteries with ST80i as they become unreliable after many charging cycles.
- **Are StressVue and ST80i studies compatible with each other?**

StressVue and ST80i studies are not compatible. StressVue and its Remote Review Station cannot read studies from the ST80i and vice versa (exception is PDFs).
- **When is the USB to Serial Port isolation adapter required?**
 - Most exercise devices require a RS232/Serial Port interface and the Tango+ BP device requires a Serial Port.
 - Many PCs today do not offer Serial Ports. Therefore, the adapter will most likely be required, ensuring proper isolation of the exercise/BP device.
- **Can the PageWriter Cardiograph PIM be used with ST80i?**
 - No, only the ST80i PIM can be used with ST80i.
- **What is the maximum distance between the AIM and PIM?**
 - The maximum distance is 3m Line of Sight (LOS).
- **How often do the HR, ST Map and Average ST complexes update on the ST80i?**

HR, ST Map and Average ST complexes update every 1 sec. on the screen.